

LOT I:
MEDICAL ITEMS AND CONSUMABLES

SN	ITEM	SPECIFICATION	UNIT OF MEASURE	QUANTITY
1	Examination Gloves	<ul style="list-style-type: none"> For medical examination Nitrile (preferable), latex, polychloroprene or PVC Powder-free Non-sterile (e.g. minimum 230 mm total length) Minimum thickness 0.05 mm. Sizes Small Available medium and Large (M, L), equally distributed. 	Box	50
2	Masks	<ul style="list-style-type: none"> Good breathability Internal and external faces should be clearly identified Efficacy passed EN 14683 Type I YY 0469 or YY/T 0969 Bacterial Filtration Efficiency (BFE): <ul style="list-style-type: none"> Must meet $\geq 95\%$ (Type I) Preferably $\geq 98\%$ (Type II / Type IIR) If bacterial droplet filtration is below 98% Or alternative equivalent standard Design and Fit-Should cover nose, mouth, and chin securely with minimal gaps., Nose clip or molding strip is recommended for better fit. 	Box	50
3	Wipes	<ul style="list-style-type: none"> Dual action cleaner and disinfectant Size – 220mm x 200mm (approx.) Proven 99.999% bacteria kill rate Very low foaming Recyclable packaging Sanitising efficacy passed & exceeds EN13727, EN1276, EN14476 Effective within 60 seconds Eliminates cross-contamination Polypropylene fabric Z-fold dispensing Plastic Sealable lid. Dual action cleaner and sanitiser Alcohol Free, Gentle on Surfaces 	pack	50
4	Sterile swabs	<ul style="list-style-type: none"> Packaged and sterilized to prevent contamination. Maximal Absorption – Engineered swab tip efficiently collects and retains contents Soft and Gentle Tip 	pcs	2000

		<ul style="list-style-type: none"> • Ergonomic Handle • Versatile Compatibility – • Ensures high reproducibility in clinical and research settings. • Packaging Options – Individually packed in sterile transport tubes, bulk options available. • Shelf Life – One month when kept frozen 		
5	10 mL Sterile Water	<ul style="list-style-type: none"> • Volume: 10 millilitres (mL) per unit • Sterility: Sterilized via autoclaving or filtration; free from microorganisms • Packaging: Single-use, sealed ampoule or vial (glass or medical-grade plastic) • Water Quality: Meets USP, EP, or applicable pharmacopeial standards (e.g., Water for Injection) • pH Range: Typically, between 5.0 – 7.0 (neutral to slightly acidic) • Endotoxin Level: Less than 0.25 EU/mL (if used for parenteral applications) • Additives: None (free from preservatives, antimicrobial agents, or buffers) • Intended Use: Suitable for dilution, reconstitution of medications, rinsing, or laboratory assays • Shelf Life: Usually 24–36 months from manufacturing date, if unopened 	pcs	200
6	Screw Cap Tubes	<ul style="list-style-type: none"> • Cap Type: Standard/ looped screw caps. • Cap Colour Options: –White. • Seal Integrity: Fitted with EPDM 'O' ring to ensure a tight, leak-proof seal, minimizing evaporation and preventing leakage. • Material: Made from high-purity virgin USP Class VI medical-grade polypropylene (PP), compliant with US FDA 21 CFR standards. • Tube Colour: Transparent for clear visibility of contents. • Capacity: 30 mL • Mouth Diameter: 25 mm (0.98 inches) • Tube Height: 90 mm (3.54 inches) • Overall Dimensions: 105 mm (Height) × 23.5 mm (Outer Diameter) / 4.13 × 0.93 inches • Packaging: Supplied in a pack of 60 screw cap tubes with lid 	pcs	400

LOT II
EQUIPMENT FOR CRIME SCENE INVESTIGATION

S/N	ITEM	DESCRIPTION/ SPECS	UNIT OF MEASURE	QUANTITY
1	Police line: Crime scene barrier tape	<ul style="list-style-type: none"> Barrier Tape black text on yellow background High quality, non-adhesive plastic tap printed with the words 'CRIME SCENE DO NOT CROSS Tape is 75mm wide x 250M long 	Roll	50
2	Tamper Evident Bags	<ul style="list-style-type: none"> Material: polyethylene (LDPE or HDPE), co-extruded for strength and durability. Seal: -Permanent adhesive or void tape closure; irreversible once opened. Thickness: -Typically ranges from 50 to 100 microns Size Options: 6"x9", 9"x12", 12"x16", etc.) to suit different contents, equally distributed. Closure Type: Single-use peel-and-seal adhesive strip or zip-lock with tamper-evident feature. Transparency: opaque Certifications: Often compliant with ISO standards or law enforcement/clinical handling guidelines. Additional Features: Write-on panels, tear-off receipts, or dual pockets (for documentation and contents separation Each bag must have a unique, pre-printed serial number and barcode Tamper-evident adhesive closure Printing Requirements: must be clearly pre-printed with: EVIDENCE" in capital letters, CHAIN OF CUSTODY" section FINAL DISPOSAL" section 	Pcs	300
3	Printed A5 Fingerprint Lift Sheets	White A5 Sheets – 80gsm - Printed with Fingerprinting Lifting Details – Supplied in plastic clip bag	pcs	50
4	Pony Hair Fingerprint brush	Pony hair fingerprint brush – Pony hair fingerprint brush 14B/1, Flat ended – Size, overall length 225mm, with pony hair length of 32mm with a width of 25mm	pcs	50

5	30g Black Powder	<ul style="list-style-type: none"> 30g Conventional Black Powder from Scene S Forensics – Supplied in screw top pot. The powder particles adhere to the latent fingerprint residue (oil, sweat, etc.). Allowing investigators to photograph and analyze them. 	bottle or jar	50
6	Fingerprint Tape	<ul style="list-style-type: none"> Latent Print Lifting Tape – 2" W x 360" L – Transparent/Tube This is a roll of adhesive tape. To prevent contamination, tape rolls can be kept in either zip lock bags or plastic tubes. Often used in the field during initial crime scene processing. 	roll	50
7	Bridge Camera	18 MP minimum - ISO 100-6400 sensitivity – Tripod mount – Min 2.5" Screen – Min 24X Zoom – Rechargeable Battery or in-built battery or comes with 4 x sets of AA rechargeable batteries (battery must fit into camera) - Charger to be 3-pin type or usb with 3-pin plug on it	pcs	1
8	Memory Card	SDHC Memory card – 16GB – Speed 10	pcs	1
9	Memory Card Reader	Memory Card Reader, BENFEI 4in1 USB 3.0 and USB-C to SD Micro SD MS CF Card Reader Adapter, 4 Cards Simultaneously Read and Write	pcs	1
10	Memory card	Plastic numbered signs (1 - 20) – Dimensions: 128 x 90 x 1.5mm – Made from 1.5mm thick vinyl material - A-frame form allows for nested storage - (+/- 10% on all sizes)	pcs	1
11	Scalpel	<ul style="list-style-type: none"> Re-usable blade Material with plastic cap: High-quality carbon or stainless steel. Length: Around 6 inches (approx. 150 mm) in total length for standard models Handle material: Textured plastic, ergonomic 	Box	20
12	Backpack	<ul style="list-style-type: none"> Responder Crime Scene backpack Capacity: 25–35 liters (sufficient for daily investigative gear) Dimensions : Approx. 45–50 cm (H) x 30–35 cm (W) x 18–22 cm (D) 	pcs	50

		<ul style="list-style-type: none"> • compartment: Large enough to hold case files, binders, and evidence kits, Secondary compartment: With multiple pockets for small equipment (e.g., flashlight, gloves, evidence markers) • Professional and low-profile collars: Black, Navy, Dark Grey • Ideal for carrying forensic and investigative gear 		
13	Plastics bag	<ul style="list-style-type: none"> ▪ Tamper-Evident Plastic Bags, Large, Plain Material: Clear or semi-transparent polyethylene (PE) or polypropylene (PP) 5cm ▪ Size (Large): Approx. 380 mm x 510 mm (<i>15 in x 20 in</i>) ▪ Thickness: Minimum 75–100 microns (heavy-duty, tear-resistant) ▪ Markings: Plain 	pcs	50
14	Scissors	160mm scissors- stainless steel with black handles	pcs	50